

Anti-Deformed Epidermal Autoregulatory Factor 1 antibody (500-565) {FITC} (STJ500749)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Deformed Epidermal Autoregulatory Factor 1 (500-565) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54006
Gene Symbol	Deaf1
Uniprot ID	DEAF1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 500-565 on mouse Deformed epidermal autoregulatory factor 1 homolog protein.
Immunogen Region	500-565
Specificity	
Immunogen Sequence	