

Anti-Disembodied antibody (STJ500763)

STJ500763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Disembodied is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IP/WB
Host/Source	Rabbit
Reactivity	D.melanogaster

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-1.2 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IHC: 1:250
	ICC: 1:250
	IF: 1:250
Formulation	CFM: 1:250
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to unique sequence on Disembodied protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	