

Anti-DNER antibody (700-737) {FITC} (STJ500779)

STJ500779

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNER (700-737) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000
	ELISA: 1:50, 000
	IP: 1:200
Formulation	IHC: 1:50-1:200
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	92737
Gene Symbol	DNER
Uniprot ID	DNER_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-737 on human Delta and Notch-like Epidermal growth factor-related receptor protein.
Immunogen Region	700-737
Specificity	
Immunogen Sequence	