

Anti-ENT1 antibody (410-460) (STJ500882)

STJ500882

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ENT1 (410-460) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	63959
Gene Symbol	Slc29a1
Uniprot ID	S29A1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 410-460 on mouse Equilibrative nitrobenzylthioinosine-sensitive nucleoside transporter ENT1 protein.
Immunogen Region	410-460
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat small intestine with Anti-ENT1 primary antibody (STJ500882) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.