

Anti-C11orf1 antibody (100-150) {Biotin} (STJ501063)

STJ501063

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C11orf1 (100-150) is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.67 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64776
Gene Symbol	CFAP68
Uniprot ID	CFA68_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human C11orf1 protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	