

Rabbit Anti-Human/Mouse/Rat Light Chain antibody {Biotin} (STJ501176)

STJ501176

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated rabbit polyclonal anti-Human/Mouse/Rat Light Chain secondary antibody. For use in most research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Rabbit anti-Goat IgG (L)
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	