

Anti-GLUT2 antibody (473-523) (STJ501206)

STJ501206

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GLUT2 (473-523) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:8, 000
	ELISA: 1:8, 000
	IP: 1:200
	IHC: 1:100
	ICC: 1:100
	IF: 1:100
Formulation	CFM: 1:100
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25139
Gene Symbol	Slc2a4
Uniprot ID	GLUT4_RAT
Immunogen	Synthetic peptide taken within amino acid region 473-523 on mouse GLUT2 protein.
Immunogen Region	473-523
Specificity	
Immunogen Sequence	