

Anti-SLC2A5 antibody {FITC} (STJ501222)

STJ501222

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC2A5 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Chimpanzee/Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:1, 000 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6518
Gene Symbol	SLC2A5
Uniprot ID	GTR5_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human GLUT 5 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	