

Anti-GOLGA4 antibody (2025-2075) (STJ501241)

STJ501241

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GOLGA4 (2025-2075) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2803
Gene Symbol	GOLGA4
Uniprot ID	GOLGA4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2025-2075 on human Golgin-245 protein.
Immunogen Region	2025-2075
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse kidney with Anti-GOLGA4 primary antibody (STJ501241) at 1:90 dilution in Buffer. Section was treated with Haematoxylin stain.