

Anti-GPR75 antibody (1-50) (STJ501262)

STJ501262

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GPR75 (1-50) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IHC/IF/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-1.25 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10936
Gene Symbol	GPR75
Uniprot ID	GPR75_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human G protein-coupled receptor 75 protein.
Immunogen Region	1-50
Specificity	Will detect extracytoplasmic domain of GPCR75 protein.
Immunogen Sequence	



Western blot with Anti-GPR75 Antibody (STJ501262) of Human RPE cells at 1:500 dilution in buffer.