

Anti-GPR75 antibody (25-75) {FITC} (STJ501268)

STJ501268

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GPR75 (25-75) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IHC/IF/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-1.25 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IP: 1:250
	IHC: 1:250
	ICC: 1:250
Formulation	IF: 1:250
	CFM: 1:250
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10936
Gene Symbol	GPR75
Uniprot ID	GPR75_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 25-75 on human G protein-coupled receptor 75 protein.
Immunogen Region	25-75
Specificity	
Immunogen Sequence	