

Anti-Hepatitis C virus Core NS3 + NS4 antibody (C-Term) (STJ501317)

STJ501317

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hepatitis C virus Core NS3 + NS4 (C-Term) is suitable for use in Dot Blot, ELISA and Western Blot research applications.
Applications	DB/ELISA/WB
Host/Source	Rabbit
Reactivity	Hepatitis C

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.65 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	C-epitope specific synthetic peptide corresponding to unique amino acid sequences on hepatitis C polyprotein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	