

Anti-HIPK2 antibody (1150-1198) {Biotin} (STJ501344)

STJ501344

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HIPK2 (1150-1198) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	28996
Gene Symbol	HIPK2
Uniprot ID	HIPK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1150-1198 on human Homeodomain-interacting protein kinase 2 protein.
Immunogen Region	1150-1198
Specificity	
Immunogen Sequence	