

## Anti-Outermembrane protein 26 H. pylori antibody (STJ501361)

STJ501361

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Outermembrane protein 26 H. pylori is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	CM/ELISA/ICC/IF/IP/IHC/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	H.pylori

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.80-1.35 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution</b>	WB: 1:500
<b>Range</b>	ELISA: 1:5, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:300 CFM: 1:200
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Heat killed and denatured Helicobacter pylori extracts.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	