

Anti-hopQ H. Pylori antibody {FITC} (STJ501366)
STJ501366**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-hopQ H. Pylori is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IP/IHC/WB
Host/Source	Rabbit
Reactivity	H.pylori

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.80-1.35 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:5, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:300 CFM: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Specificity	
Immunogen Sequence	