

Anti-OPTN antibody {Biotin} (STJ501368)

STJ501368

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OPTN is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IP/IHC/WB
Host/Source	Rabbit
Reactivity	H.pylori

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.80-1.35 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:500
Range	ELISA: 1:5, 000
	IP: 1:150
	IHC: 1:200
	ICC: 1:200
	IF: 1:300
Formulation	CFM: 1:200
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10133
Gene Symbol	OPTN
Uniprot ID	OPTN_HUMAN
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	