

Anti-IRF8 antibody (350-400) (STJ501500)

STJ501500

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IRF8 (350-400) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	15900
Gene Symbol	Irf8
Uniprot ID	IRF8_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 350-400 on mouse Interferon Regulatory Factor 8 protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	



Western blot with Anti-IRF8 Antibody (STJ501500) with
A) RAW cells, B) Mouse liver at 1:500 dilution in buffer.