

Anti-LOXL1 antibody (570-607) {FITC} (STJ501635)

STJ501635

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LOXL1 (570-607) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500-1:1, 000
Range	ELISA: 1:20, 000-1:50, 000 IHC: 1:50 IF: 1:50
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16949
Gene Symbol	Loxl1
Uniprot ID	LOXL1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 570-607 on rat LOXL1 protein.
Immunogen Region	570-607
Specificity	
Immunogen Sequence	