

**Anti-MANF antibody (80-120) {Biotin} (STJ501655)**

STJ501655

**GENERAL INFORMATION**

|                          |                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                            |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-MANF (80-120) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | ELISA/IHC/IP/WB                                                                                                                                               |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                        |
| <b>Reactivity</b>        | Human/Mouse/Rat                                                                                                                                               |

**PRODUCT PROPERTIES**

|                       |                                                                                                                                                                      |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                                                                           |
| <b>Clone ID</b>       |                                                                                                                                                                      |
| <b>Concentration</b>  | 0.5 µg/µl                                                                                                                                                            |
| <b>Conjugation</b>    | Biotin                                                                                                                                                               |
| <b>Purification</b>   | Affinity Purified                                                                                                                                                    |
| <b>Dilution Range</b> | WB: 1:500<br>ELISA: 1:50, 000<br>IP: 1:200<br>IHC: 1:50-1:150                                                                                                        |
| <b>Formulation</b>    | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| <b>Isotype</b>        | IgG                                                                                                                                                                  |
| <b>Storage</b>        | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |
| <b>Instruction</b>    |                                                                                                                                                                      |

**TARGET INFORMATION**

|                           |                                                                                                                       |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 7873                                                                                                                  |
| <b>Gene Symbol</b>        | MANF                                                                                                                  |
| <b>Uniprot ID</b>         | MANF_HUMAN                                                                                                            |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 80-120 on human Mesencephalic astrocyte-derived neurotrophic factor. |
| <b>Immunogen Region</b>   | 80-120                                                                                                                |
| <b>Specificity</b>        |                                                                                                                       |
| <b>Immunogen Sequence</b> |                                                                                                                       |