

Anti-MYLK antibody (1870-1907) {FITC} (STJ501836)

STJ501836

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-MYLK (1870-1907) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgY
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4638
Gene Symbol	MYLK
Uniprot ID	MYLK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1870-1907 on human Myosin light chain kinase, smooth muscle isoform 1 protein.
Immunogen Region	1870-1907
Specificity	
Immunogen Sequence	