

Anti-NPAS4 antibody (STJ501951)

STJ501951

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPAS4 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1;10, 000 ELISA: 1;10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	266743
Gene Symbol	NPAS4
Uniprot ID	NPAS4_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human Npas4 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-NPAS4 Antibody (STJ501951) with peptide on at 1:500 dilution in buffer.