

Anti-OLIG2 antibody (136-156) {Biotin} (STJ502006)

STJ502006

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OLIG2 (136-156) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10215
Gene Symbol	OLIG2
Uniprot ID	OLIG2_HUMAN
Immunogen	Synthetic peptide corresponding to Human Olig2 amino acids 136-156. The selected peptide was post-synthetically modified to achieve highest antigenicity before coupling to a carrier protein using hetero-bifunctional cross linker for immunogen preparation.
Immunogen Region	136-156
Specificity	
Immunogen Sequence	