

Anti-Phospho-BTK antibody (544-561) {FITC} (STJ502178)

STJ502178

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-BTK (544-561) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Human/Monkey/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	695
Gene Symbol	BTK
Uniprot ID	BTK_HUMAN
Immunogen	Synthetic peptide corresponding to Human BTK amino acids 544-561 (phospho Y551). Sequence RYV LDD EY (p) T SSV GSK FPV
Immunogen Region	544-561
Specificity	Does not cross react to any other cellular protein.
Immunogen Sequence	