

Anti-PDE3A antibody (STJ502263)

STJ502263

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-PDE3A is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Goat
Reactivity	Human/Monkey/Mouse/Rat/Xenopus

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.72 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5139
Gene Symbol	PDE3A
Uniprot ID	PDE3A_HUMAN
Immunogen	Synthetic peptide selective only for PDE3A variants
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-PDE3A Antibody (STJ502263) with PDE3A peptide. Antibody at 1:250 dilution in buffer.