

Anti-Phospho-FGFR1 antibody (640-690) (STJ502423)

STJ502423

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Phospho-FGFR1 (640-690) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA/IP/WB |
| Host/Source | Rabbit |
| Reactivity | Human/Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.75 µg/µl |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:10, 000 IP: 1:150 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 2260 |
| Gene Symbol | FGFR1 |
| Uniprot ID | FGFR1_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 640-690 on human basic fibroblast growth factor receptor protein 1. |
| Immunogen Region | 640-690 |
| Specificity | |
| Immunogen Sequence | |