

Anti-PGAP3 antibody (1-50) {FITC} (STJ502442)

STJ502442

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PGAP3 (1-50) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB/ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1, 000 DB: 1:50, 000 ELISA: 1:50, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93210
Gene Symbol	PGAP3
Uniprot ID	PGAP3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human PGAP3 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	