

**Anti-PGAP3 antibody (10-60) {Biotin} (STJ502444)**

STJ502444

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PGAP3 (10-60) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	DB/ELISA/IHC/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Monkey

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.50-0.75 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:1, 000 DB: 1:50, 000 ELISA: 1:50, 000 IP: 1:150 IHC: 1:100
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">93210</a>
<b>Gene Symbol</b>	<a href="#">PGAP3</a>
<b>Uniprot ID</b>	<a href="#">PGAP3_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 10-60 on human PGAP3 protein.
<b>Immunogen Region</b>	10-60
<b>Specificity</b>	
<b>Immunogen Sequence</b>	