

Anti-Phospho-JAK1 antibody (1000-1050) {FITC} (STJ502509)
STJ502509**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-JAK1 (1000-1050) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3716
Gene Symbol	JAK1
Uniprot ID	JAK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1050 on human JAK1 protein. Phosphorylated 1034 & 1035 tyrosine residues.
Immunogen Region	1000-1050
Specificity	
Immunogen Sequence	