

Anti-Phospho-MERTK antibody (746-757) {Biotin} (STJ502559)

STJ502559

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-MERTK (746-757) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/FC/WB
Host/Source	Rabbit
Reactivity	Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.63 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:750 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65037
Gene Symbol	Mertk
Uniprot ID	MERTK_RAT
Immunogen	Synthetic peptide corresponding to amino acids 746-757 of Human MERTK.
Immunogen Region	746-757
Specificity	
Immunogen Sequence	