

Anti-PODXL antibody (510-560) {FITC} (STJ502578)

STJ502578

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PODXL (510-560) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB/ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5420
Gene Symbol	PODXL
Uniprot ID	PODXL_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 510-560 on human PODXL protein.
Immunogen Region	510-560
Specificity	This antibody will react with both PODXL isoform 1 and PODXL isoform 2.
Immunogen Sequence	