

Anti-POP4 antibody (170-220) {FITC} (STJ502590)

STJ502590

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-POP4 (170-220) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Bovine/Cat/Human/Monkey/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10775
Gene Symbol	POP4
Uniprot ID	RPP29_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 170-220 on human Ribonuclease P protein subunit p29 protein.
Immunogen Region	170-220
Specificity	
Immunogen Sequence	