

Anti-Phospho-PDE3A antibody (200-300) {FITC} (STJ502605)

STJ502605

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PDE3A (200-300) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.0-1.25 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54611
Gene Symbol	Pde3a
Uniprot ID	PDE3A_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-300 on mouse PDE3A protein. Phosphorylated Serine 290 & 292 residues.
Immunogen Region	200-300
Specificity	This antibody will detect PDE3A protein that has been phosphorylated at its putative PKB/Akt site.
Immunogen Sequence	