

Anti-DNA polymerase eta antibody (200-300) (STJ502747)

STJ502747

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNA polymerase eta (200-300) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	<i>S.cerevisiae</i>

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.70-0.98 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 200-300 on Chain A, Yeast DNA polymerase Eta.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	