

Anti-RIP140 antibody (50-100) {FITC} (STJ502793)

STJ502793

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIP140 (50-100) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IHC/IF/IP/WB
Host/Source	Rabbit
Reactivity	Chicken/Human/Monkey/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.60 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268903
Gene Symbol	Nrip1
Uniprot ID	NRIP1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse receptor interacting protein 140.
Immunogen Region	50-100
Specificity	Recognizes RIP140 at the N terminal region containing the HDAC binding domain.
Immunogen Sequence	