

Anti-RIPK1 antibody (620-670) {FITC} (STJ502800)

STJ502800

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIPK1 (620-670) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8737
Gene Symbol	RIPK1
Uniprot ID	RIPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 620-670 on human Receptor-interacting Serine/Threonine-protein kinase 1 protein.
Immunogen Region	620-670
Specificity	
Immunogen Sequence	