

Anti-RIPK3 antibody (468-518) (STJ502805)

STJ502805

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIPK3 (468-518) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Bovine/Human/Monkey/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66-0.68 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11035
Gene Symbol	RIPK3
Uniprot ID	RIPK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 468-518 on human Receptor-interacting Serine/Threonine-protein kinase 3 protein.
Immunogen Region	468-518
Specificity	
Immunogen Sequence	