

Anti-SERHL2 antibody {FITC} (STJ502946)

STJ502946

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SERHL2 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.40-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000-1:100, 000 ELISA: 1:10, 000-1:100, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253190
Gene Symbol	SERHL2
Uniprot ID	SEHL2_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human serine hydrolase 1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	