

Anti-C8orf44-SGK3 antibody (440-490) (STJ502984)

STJ502984

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-C8orf44-SGK3 (440-490) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/IC/IF/IHC/IP/WB
Host/Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100533105/23678
Gene Symbol	SGK3
Uniprot ID	SGK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 440-490 on human Sgk3 protein.
Immunogen Region	440-490
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse small intestine with Anti-C8orf44-SGK3 primary antibody (STJ502984) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081