

Anti-SPACA3 antibody (80-130) {Biotin} (STJ503102)

STJ503102

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPACA3 (80-130) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124912
Gene Symbol	SPACA3
Uniprot ID	SACA3_HUMAN
Immunogen	Synthetic peptides taken within amino acid region 80-130 on human sperm acrosome membrane-associated protein 3 isoform 2.
Immunogen Region	80-130
Specificity	This antibody will label all isoforms of sperm acrosome membrane-associated protein 3.
Immunogen Sequence	