

Anti-DNTT antibody (459-509) {FITC} (STJ503238)

STJ503238

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNTT (459-509) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Chicken/Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:1, 000
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IP: 1:250
	IHC: 1:250
	ICC: 1:250
	IF: 1:250
Formulation	CFM: 1:250
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1791
Gene Symbol	DNTT
Uniprot ID	TDT_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 459-509 on human Terminal deoxynucleotidyltransferase protein.
Immunogen Region	459-509
Specificity	
Immunogen Sequence	