

Anti-TIA1 antibody (414-464) {FITC} (STJ503256)

STJ503256

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TIA1 (414-464) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IHC: 1:50
	ICC: 1:50
	IF: 1:50
Formulation	CFM: 1:50
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7072
Gene Symbol	TIA1
Uniprot ID	TIA1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 414-464 on human TIA1 protein.
Immunogen Region	414-464
Specificity	
Immunogen Sequence	