

## Anti-TLR11 antibody (430-460) (STJ503275)

STJ503275

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TLR11 (430-460) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Chimpanzee/Human/Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.56-0.72 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	239081
<b>Gene Symbol</b>	Tlr11
<b>Uniprot ID</b>	TLR11_MOUSE
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 430-460 on human TLR11 protein.
<b>Immunogen Region</b>	430-460
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot with Anti-TLR11 Antibody (STJ503275) with TLR11 peptide. Antibody at 1:500 dilution in buffer.