

Anti-Timeless antibody (1178-1197) {Biotin} (STJ503334)

STJ503334

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Timeless (1178-1197) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.98-1.26 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	21853
Gene Symbol	Timeless
Uniprot ID	TIM_MOUSE
Immunogen	Synthetic peptide: SGTPRVHRKKRFQIEDED corresponding to amino acids 1178-1197 on Mouse Timeless protein.
Immunogen Region	1178-1197
Specificity	
Immunogen Sequence	