

Anti-TPCN2 antibody (205-255) {FITC} (STJ503371)

STJ503371

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-TPCN2 (205-255) is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA/WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.57 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 219931 |
| Gene Symbol | TPCN2 |
| Uniprot ID | TPC2_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 205-255 on human TPC2 protein. |
| Immunogen Region | 205-255 |
| Specificity | |
| Immunogen Sequence | |