

Anti-TTC27 antibody (100-150) {FITC} (STJ503435)

STJ503435

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-TTC27 (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA/IP/WB |
| Host/Source | Rabbit |
| Reactivity | Zebra Fish |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.50 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:200 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 74196 |
| Gene Symbol | Ttc27 |
| Uniprot ID | TTC27_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 100-150 on Tetratricopeptide repeat protein 27 Danio rerio protein. |
| Immunogen Region | 100-150 |
| Specificity | |
| Immunogen Sequence | |