

Anti-ULK1 antibody (1000-1050) {Biotin} (STJ503462)

STJ503462

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ULK1 (1000-1050) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Bovine/Human/Monkey/Mouse/Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:250 ELISA: 1:10, 000 IP: 1:200 IF: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8408
Gene Symbol	ULK1
Uniprot ID	ULK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1050 on mouse ULK1 protein.
Immunogen Region	1000-1050
Specificity	
Immunogen Sequence	