

Anti-USH2A antibody (1495-1545) (STJ503475)

STJ503475

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-USH2A (1495-1545) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:4, 000
	ELISA: 1:4, 000
	IP: 1:200
	IHC: 1:250
	ICC: 1:250
	IF: 1:250
Formulation	CFM: 1:250
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7399
Gene Symbol	USH2A
Uniprot ID	USH2A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1495-1545 on human Usherin isoform A protein.
Immunogen Region	1495-1545
Specificity	This antibody will detect both Usherin isoform A and Usherin isoform B.
Immunogen Sequence	