

Anti-VKORC1 antibody (10-60) {FITC} (STJ503490)

STJ503490

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VKORC1 (10-60) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.70 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:50, 000
	ELISA: 1:50, 000
	IP: 1:200
	IHC: 1:250
	ICC: 1:250
Formulation	IF: 1:250
	CFM: 1:250
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	309004
Gene Symbol	Vkorc1
Uniprot ID	VKOR1_RAT
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human VKORC1 protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	