

Anti-WWOX antibody (250-300) {Biotin} (STJ503544)

STJ503544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WWOX (250-300) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Monkey/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000-1:100, 000 ELISA: 1:10, 000-1:100, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51741
Gene Symbol	WWOX
Uniprot ID	WWOX_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human WWOXdelta5-8 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	