

Anti-Xflop antibody (200-250) {FITC} (STJ503548)

STJ503548

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Xflop (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Xenopus

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-1.20 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 200-250 on Xenopus laevis orphan G protein-coupled receptor Xflop protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	